

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

### 1 Identification of the substance/mixture and of the company/undertaking

- **Product identifier**
- **Trade name:** PrimaDry™
- **Article number:** SDS 65-001.12R01, 10955, 10954, 10993, 716, 716-JP, 725-JP, 7663
- **Relevant identified uses of the substance or mixture and uses advised against**  
Professional Dental Priming Agent
- **Application of the substance / the mixture** Professional Dental Priming Agent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ultradent Products Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com  
(800) 552-5512
- **EC Responsible Person**  
Ultradent Products GmbH  
Am Westhover Berg 30  
51149 Cologne Germany  
Email: infoDE@ultradent.com  
Office Phone: +49(0)2203-35-92-0
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**  
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

### 2 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



Flam. Liq. 2 H225 Highly flammable liquid and vapour.



Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2 H319 Causes serious eye irritation.

- **Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Void
- **Hazard pictograms** GHS02, GHS07
- **Signal word** Danger
- **Hazard statements**  
H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.

(Contd. on page 2)

# Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 1)

*H319 Causes serious eye irritation.***· Precautionary statements***P101 If medical advice is needed, have product container or label at hand.**P102 Keep out of reach of children.**P103 Read carefully and follow all instructions.**P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.**P241 Use explosion-proof [electrical/ventilating/lighting] equipment.**P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].**P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.**P403+P235 Store in a well-ventilated place. Keep cool.**P501 Dispose of contents/container in accordance with local/regional/national/international regulations.*

### 3 Composition/information on ingredients

**· Mixtures****· Description:** Mixture of substances listed below with nonhazardous additions.**· Dangerous components:**

CAS: 64-17-5 EINECS: 200-578-6	Ethyl Alcohol	 Flam. Liq. 2, H225	<100%
	Trade Secret	 Skin Corr. 1A, H314; Eye Dam. 1, H318	>1-<3%

**· Additional information:** For the wording of the listed hazard phrases refer to section 16.

### 4 First aid measures

**· Description of first aid measures****· General information:**

*Take proper precautions to ensure your own health and safety before attempting rescue and providing first aid. Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area. Immediately remove any clothing soiled by the product.*

**· After inhalation:***Seek medical treatment in case of complaints.**Give oxygen or artificial respiration as needed.**Supply fresh air.**In case of unconsciousness place patient stably in side position for transportation.***· After skin contact:***If skin irritation continues, consult a doctor.**Launder clothing before reuse.**Immediately wash with water and soap and rinse thoroughly.***· After eye contact:***Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.***· After swallowing:***Do NOT induce vomiting.**If vomiting does occur, have victim lean forward to prevent aspiration.**Rinse mouth with water.**Seek medical treatment.**Never give anything by mouth to an unconscious person.*

(Contd. on page 3)

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 2)

- **Most important symptoms and effects, both acute and delayed**  
Symptoms vary with alcohol level of the blood. Mild alcohol intoxication occurs at blood levels between 0.05-0.15%. Approximately 25% of individuals show signs of intoxication at these levels. Above 0.15% the person is definitely under the influence of ethanol; 50-95% of individuals are clinically intoxicated at these levels. Severe poisoning occurs when the blood ethanol level is 0.3- 0.5%. Above 0.5% the individual will be comatose and death can occur. The unabsorbed ethanol should be removed by gastric lavage after intubating the patient to prevent aspiration. Avoid the use of depressant drugs or the excessive administration of fluids.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
Water fog  
Alcohol resistant foam  
Water spray  
Carbon dioxide  
Dry Chemical
- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **Special hazards arising from the substance or mixture**  
Carbon monoxide (CO)  
May produce a floating fire hazard.  
Static ignition hazard can result from handling and use.  
Vapors may travel to source of ignition and flash back.  
Vapors may settle in low or confined spaces.  
  
Alcohols burn with a pale blue flame which may be extremely hard to see under normal lighting conditions. Personnel may only be able to feel the heat of the fire without seeing flames. Extreme caution must be exercised in fighting alcohol fires. Fight fire from maximum distance or use unmanned hose holders or monitor nozzles. Cool containers with flooding quantities of water until well after fire is out. Withdraw immediately in case of rising sound from venting safety devices or discoloration of tank. Always stay away from tanks engulfed in fire.
- **Advice for firefighters:**
- **Protective equipment:**  
Wear self-contained respiratory protective device.  
Wear fully protective suit.
- **Additional information** Cool endangered receptacles with water spray.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**  
Do not inhale vapors, mist or gas. Ensure adequate ventilation. Remove all sources of ignition. Evacuate personnel to safe areas. Beware of vapors accumulating to form explosive concentrations. Vapors can accumulate in low areas.  
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:**  
Stop leak. Contain spill if possible and safe to do so.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Highly flammable liquid. Eliminate all sources of ignition. All equipment used when handling this product must be grounded. A vapor suppressing foam may be used to reduce vapors. Do not touch or walk through spilled material. Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations.

(Contd. on page 4)

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

Trade name: PrimaDry™

(Contd. of page 3)

- Use clean non-sparking tools to collect absorbed material.
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- Dispose contaminated material as waste according to section 13.
- Ensure adequate ventilation.
- **Reference to other sections**
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

### 7 Handling and storage

- **Precautions for safe handling:**
- Avoid contact with eyes, skin, and clothing.
- Do not inhale vapor or mist.
- **Information about fire - and explosion protection:**
- Keep ignition sources away - Do not smoke.
- Protect against electrostatic charges.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**
- Metal containers involved in the transfer of this material should be grounded and bonded.
- Store in a cool location.
- **Information about storage in one common storage facility:** Store away from flammable substances.
- **Further information about storage conditions:**
- Containers which are opened must be carefully resealed and kept upright to prevent leakage.
- Consult local fire codes for additional storage information.
- See product labelling.
- Keep container tightly sealed.
- Store in cool, dry conditions in well - sealed receptacles.
- **Specific end use(s)** Professional Dental Priming Agent

### 8 Exposure controls/personal protection

- **Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**
- 64-17-5 Ethyl Alcohol**
- WEL Long-term value: 1920 mg/m<sup>3</sup>, 1000 ppm
- **Additional information:** The lists valid during the making were used as basis.
- **Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
- Observe good industrial hygiene practices.
- Ensure that washing facilities are available at the work place.
- Electrical equipment should be grounded and confirm to applicable electrical code.
- When using do not smoke.
- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes and skin.

(Contd. on page 5)

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 4)

· **Respiratory protection:**

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

· **Hand protection**

Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product.

Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye/face protection**

Maintain eye wash fountain and quick-drench facilities in work area.

Use equipment approved by appropriate government standards, such as NIOSH (US) or EN166 (EU)

Face protection

Use chemical safety goggles and/or a full face shield where splashing is possible.



Tightly sealed goggles

· **Body protection:**

Choose body protection according to the amount and concentration of the dangerous substance at the work place.

### 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Physical state**

Fluid

· **Colour:**

Yellow tint

· **Odour:**

Alcohol-like

· **Odour threshold:**

Not determined.

· **Melting point/freezing point:**

Undetermined.

· **Boiling point or initial boiling point and boiling range**

78 °C

· **Flammability**

Highly flammable.

· **Lower and upper explosion limit**

· **Lower:**

3.3 Vol %

· **Upper:**

19 Vol %

(Contd. on page 6)

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 5)

· <b>Flash point:</b>	13 °C
· <b>Auto-ignition temperature:</b>	363 °C
· <b>Decomposition temperature:</b>	Not determined.
· <b>pH</b>	Not determined.
· <b>Viscosity:</b>	
· <b>Kinematic viscosity</b>	Not determined.
· <b>Dynamic:</b>	Not determined.
· <b>Solubility</b>	
· <b>water at 20 °C:</b>	1,000 g/l
· <b>Partition coefficient n-octanol/water (log value)</b>	Not determined.
· <b>Vapour pressure:</b>	Not determined.
· <b>Density and/or relative density</b>	
· <b>Density at 20 °C:</b>	0.79 g/cm <sup>3</sup>
· <b>Relative density</b>	Not determined.
· <b>Vapour density at 20 °C</b>	1.6 g/cm <sup>3</sup>

· <b>Other information</b>	
· <b>Appearance:</b>	
· <b>Form:</b>	Liquid
· <b>Important information on protection of health and environment, and on safety.</b>	
· <b>Ignition temperature:</b>	Product is not selfigniting.
· <b>Explosive properties:</b>	Product is not explosive. However, formation of explosive air/vapour mixtures are possible.
· <b>Change in condition</b>	
· <b>Evaporation rate</b>	Not determined.

· <b>Information with regard to physical hazard classes</b>	
· <b>Explosives</b>	Void
· <b>Flammable gases</b>	Void
· <b>Aerosols</b>	Void
· <b>Oxidising gases</b>	Void
· <b>Gases under pressure</b>	Void
· <b>Flammable liquids</b>	Highly flammable liquid and vapour.
· <b>Flammable solids</b>	Void
· <b>Self-reactive substances and mixtures</b>	Void
· <b>Pyrophoric liquids</b>	Void
· <b>Pyrophoric solids</b>	Void
· <b>Self-heating substances and mixtures</b>	Void
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	Void
· <b>Oxidising liquids</b>	Void
· <b>Oxidising solids</b>	Void
· <b>Organic peroxides</b>	Void
· <b>Corrosive to metals</b>	Void
· <b>Desensitised explosives</b>	Void

### 10 Stability and reactivity

- **Reactivity Stable**
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.

(Contd. on page 7)

# Safety data sheet

## according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

Trade name: PrimaDry™

(Contd. of page 6)

- **Possibility of hazardous reactions:**  
Vapors may form explosive mixture with air.  
Flammable.
- **Conditions to avoid:**  
Direct sunlight  
Extreme temperature  
Flames  
Sparks  
Heat
- **Incompatible materials:**  
Alkali metals  
Strong Inorganic Acids  
Peroxides  
Oxidizing Agents  
Ammonia
- **Hazardous decomposition products:** Carbon monoxide and carbon dioxide

### 11 Toxicological information

- **Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.

- **LD/LC50 values relevant for classification:**

#### 64-17-5 Ethyl Alcohol

Oral	LD50	5,600 mg/kg (guinea pig)
		3,400 mg/kg (mouse)
		7,060 mg/kg (rat)
Inhalative	LC50 Fish	>10,000 mg/l (Fish)
	LC50/4 h	39 mg/l (mouse) 20,000 mg/l (rat)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Causes skin irritation.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Information on other hazards**

- **Endocrine disrupting properties**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**

- **Aquatic toxicity:**

#### 64-17-5 Ethyl Alcohol

Algae Toxicity | 1,000 mg/l (Algae)

- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

(Contd. on page 8)

## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

Trade name: PrimaDry™


(Contd. of page 7)

- **Endocrine disrupting properties** The product does not contain substances with endocrine disrupting properties.
- **Other adverse effects**
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |   |   |
|---|---|
| <b>UN number or ID number</b>   | UN1987  |
| <b>ADR, IMDG, IATA</b>  |   |
| <b>UN proper shipping name</b>  | 1987 ALCOHOLS, N.O.S. (Ethyl Alcohol)   |
| <b>ADR</b>  | ALCOHOLS, N.O.S. (Ethyl Alcohol)  |
| <b>IMDG, IATA</b>   |   |
| <b>Transport hazard class(es)</b>   |   |
| <b>ADR, IMDG, IATA</b>  |   |
|  |   |
| <b>Class</b>  | 3 Flammable liquids.  |
| <b>Label</b>  | 3   |
| <b>Packing group</b>  | II  |
| <b>ADR, IMDG, IATA</b>  |   |
| <b>Environmental hazards:</b>   | Not applicable.   |
| <b>Special precautions for user</b>   | Warning: Flammable liquids.   |
| <b>Hazard identification number (Kemler code):</b>                                  | 33  |
| <b>EMS Number:</b>  | F-E,S-D   |
| <b>Stowage Category</b>   | B   |
| <b>Maritime transport in bulk according to IMO instruments</b>                      | Not applicable.   |
| <b>Transport/Additional information:</b>  |   |
| <b>ADR</b>  |   |
| <b>Limited quantities (LQ)</b>  | 1L  |
| <b>Excepted quantities (EQ)</b>   | Code: E2<br>Maximum net quantity per inner packaging: 30 ml<br>Maximum net quantity per outer packaging: 500 ml |
| <b>Transport category</b>   | 2   |

(Contd. on page 9)



## Safety data sheet according to UK REACH

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 8)

· <b>Tunnel restriction code</b>	D/E
· <b>IMDG</b>	
· <b>Limited quantities (LQ)</b>	1L
· <b>Excepted quantities (EQ)</b>	Code: E2 Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 500 ml
· <b>UN "Model Regulation":</b>	UN 1987 ALCOHOLS, N.O.S. (ETHYL ALCOHOL), 3, II

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **Poisons Act**

· **Regulated explosives precursors**

None of the ingredients is listed.

· **Regulated poisons**

None of the ingredients is listed.

· **Reportable explosives precursors**

None of the ingredients is listed.

· **Reportable poisons**

None of the ingredients is listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

· **Seveso category P5c** FLAMMABLE LIQUIDS

· **Qualifying quantity (tonnes) for the application of lower-tier requirements** 5,000 t

· **Qualifying quantity (tonnes) for the application of upper-tier requirements** 50,000 t

· **Chemical safety assessment:**

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases from Section 3**

H225 Highly flammable liquid and vapour.

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

(Contd. on page 10)

**Safety data sheet**  
**according to UK REACH**

Printing date 12.09.2025

Version number 1

Revision: 12.09.2025

**Trade name: PrimaDry™**

(Contd. of page 9)

*CAS: Chemical Abstracts Service (division of the American Chemical Society)**LC50: Lethal concentration, 50 percent**LD50: Lethal dose, 50 percent**PBT: Persistent, Bioaccumulative and Toxic**vPvB: very Persistent and very Bioaccumulative**NIOSH: National Institute for Occupational Safety**ATE: Acute toxicity estimate values**Flam. Liq. 2: Flammable liquids – Category 2**Skin Corr. 1A: Skin corrosion/irritation – Category 1A**Skin Irrit. 2: Skin corrosion/irritation – Category 2**Eye Dam. 1: Serious eye damage/eye irritation – Category 1**Eye Irrit. 2: Serious eye damage/eye irritation – Category 2**· \* **Data compared to the previous version altered.***

GB